

Healthcare
furniture,
enhancing
lives.

Caremed **Alrick.**

2400 Series
Bariatric Bed




The Alrick 2400 Series Bariatric Tilt Bed provides the ultimate comfort for bariatric patients with its safe working load reaching 400kg. With its features such as 150mm Tente Twin Wheel Central Lock Castors it allows the large bed to be easily and safely manouvered into position. This bed is now available in King Single, Extra Wide, Double & Queen.


Configuration Possibilities




Accessories Available




Standard Functions




Electric Hi-Low



Electric Backrest



Electric Forward Trendelenburg



Electric Reverse Trendelenburg

Standard features also include:

- 150mm twin wheel central locking castors
- Battery backup
- Underbed lighting
- Backlit handset
- CPR Emergency release
- APO Grey top frame powdercoat
- Ironstone base frame & endboard frame powdercoat
- Fini Blonde laminate head and footboards

Specifications

NET Bed weight	King Single	120kg
Low Position	390mm	
High Position	830mm	
Underbed clearance	160mm	
Angles (°)	Back	65
	Knee	44
Warranty	Frame	10 yr
	Electric	5 yr
	Castors	1 yr
	Handset	1 yr
Optional features	ACO Attendant Control	
Product codes	2400KS-T	King Single
	2400EW-T	Extra Wide
	2400D-T	Double
	2400Q-T	Queen
	2400K-T	King Double

	Overall dimensions	Sleeping surface dimensions	SWL
King Single	2170 x 1090mm	2000 x 1070mm	400kg
Extra Wide	2170 x 1220mm	2000 x 1200mm	400kg
Double	2170 x 1390mm	2000 x 1370mm	350kg
Queen	2170 x 1520mm	2000 x 1500mm	350kg
King Double	2170 x 1820mm	2000 x 1800mm	350kg

** LOAD CAPACITY HAS BEEN TESTED IN ACCORDANCE WITH AS/NZS 3200.2.38:2007 **