

Healthcare
furniture,
enhancing
lives.

Caremed Alrick.

5000 Series Homecare Bed

Community Home Care Bed



The Alrick 5000 Pull-A-Part Transportable Series is targeted at easy transportation and storage. Ideal for short term care, the 5000 Series can be assembled and dismantled with ease through the incorporation of a simple slot system which is secured with a pin, and folds compactly for storage. Consequently the 5000 Series is an ideal homecare bed, or as a short-term care alternative.

Specifications

Height	Minimum	295mm
	Maximum	700mm
Overall Width	P5500	930mm
	P5500W	1070mm
Overall Length	2190mm	
Sleeping Surface Width	P5500	900mm
	P5500W	1050mm
Sleeping Surface Length	2000mm	
SWL	200kg	
Product Codes	Standard Single	P5500
	King Single	P5500W
Warranty	Frame	10 yr
	Electric	5 yr

Features



Pin Slot System



Transport Mode (Folded)



Handset



P5500 Standard Single
900mm Wide



P5500W King Single
1070mm Wide

Electric Hi/Low

Electric Backrest

Electric Kneebend

Manual Vascular Support

Electric Forward & Reverse
Trendelenburg


100mm Braking Castors



Australian
Made



10 year frame
warranty



5 year Linak
warranty

Options



ABB Battery Backup

Accessories




360HR Horizontal Side Rails



380VR Vertical Side Rails



400HR Deluxe Rails



ASHP/H Self Help Pole



Standezy



CLP400 Clamp On
Support Rail 400mm



COT PRO Padded
Side Rail Protectors



AIVP I.V. Pole



C.EXT5000 200mm Clip On Base Extension
C.EXT5000W 200mm Clip On Base Extension

VIC
2/77 Ricky Way
Epping VIC 3076

NSW
104 Loughnan Street
Lake Cargelligo NSW 2672

UK/EU
Port Stanley House
Falkland Close Coventry CV4 8AU

T 1800 282 001
sales@caremed.healthcare
www.caremed.healthcare