

Healthcare
furniture,
enhancing
lives.

Caremed Alrick.

EN7 Floorline Bed

It looks like it's *home*,
works like it's *in a hospital*
and lasts more than a *decade*



Need Quick Dispatch?
Contact our Team on 1800 282 001

Options

Laminate

Slender Arch Head and Foot Board



Colour: Natural Oak

Timber

Eltham V2 Head and Foot Board



Colour: Burnt Walnut - Beech - Saddle

Upholstered

Jersey Head Board and Plain Rectangle Foot Board



Colour: Marble - Ash - Ocean

Features

- Electric Hi/Low
- Electric knee bend
- Electric forward & reverse Trendelenburg
- Auto regression backrest
- 75mm Castors
- Standard head and footboards

Optional Features :

- Dual Braking
- Underbed Light
- Head and Foot Board Options
- Battery Backup
- 300kg SWL
Upgrade to 300kg SWL by changing the handset. Please note, this needs to be purchased separately.

Specifications

Overall Dimensions	Single King Single	2080 x 920mm 2080 x 1090mm
SWL		200kg
NET Bed Weight		110kg
Height	Low Position High Position	100mm 725mm

10 year frame warranty

5 year electric warranty

**Caremed
Alick.**

* Handset & Castors - 1 year

Configuration Possibilities



Code	Item Description	Colours Available
FB-EN7S0000E	EN7 Single with Eltham Headboard and Footer	Burnt Walnut, Saddle & Beech
FB-EN7S0000JU	EN7 Single with Jersey Upholstered Headboard and Plain Upholstered Footer	Ocean, Ash & Marble
FB-EN7S0000SAL	EN7 Single with Laminate Slender Arch Headboard and Footer	Natural Oak
FB-EN7KS0000E	EN7 King Single with Eltham Headboard and Footer	Burnt Walnut, Saddle & Beech
FB-EN7KS0000JU	EN7 King Single with Jersey Upholstered Headboard and Plain Upholstered Footer	Ocean, Ash & Marble
FB-EN7KS0000SAL	EN7 Single with Laminate Slender Arch Headboard and Footer	Natural Oak

VIC
2/77 Ricky Way
Epping VIC 3076

NSW
104 Loughnan Street
Lake Cargelligo NSW 2672

NSW
Level 6, 10 Herb Elliot Avenue,
Sydney Olympic Park NSW, 2127

T 1800 282 001
sales@caremed.healthcare
www.caremed.healthcare