

Air Comfort Chair

- Classic
- Slimline
- Bariatric



A classic pressure care chair featuring a seat air cushion, designed for use in high care environments. Available in Classic, Slimline and Bariatric sizes, to suit most users.

Product Specifications



Air Chair Classic Air Chair Slimline Air Chair Bariatric

Classic

Order code: ACC

Safe working load: 150kg
Overall height: 1300mm
Overall width: 710mm
Overall depth: 1250mm
Seat level height: 700mm
Seat width: 510mm
Seat depth: 560mm
Seat back adjustment angle: 97-136 °
Seat tilt adjustment angle: 0-17 °
Leg rest adjustment angle: 0-90 °

Slimline

Order code: ACS

As above, except:
Overall width: 650mm
Seat width: 450mm

Bariatric

Order code: ACB

Safe working load: 200kg
Overall height: 1300mm
Overall width: 800mm
Overall depth: 1250mm
Seat level height: 700mm
Seat width: 600mm
Seat depth: 560mm
Angle adjustments as for Classic

Unit weight: Classic/Slimline 50kg, Bariatric 60kg

7

The Air Chair is specifically designed for use in high care environments where patients can be at risk of pressure sores and other pressure management issues. It uses an air cushion in the seat and backrest area, along with advanced foam system and four-way stretch material throughout, to offer excellent pressure relief for those at risk.

To cater for a range of users and requirements, the chair is highly adjustable, including arm rests and head wings that can be lowered below seat level for side transfers, as well as back rest, leg rest and tilt-in-space adjustment.

- Specialised pressure relief and positioning support for elderly and special needs users
- Cushion air cell system in seat and backrest is pressure adjustable for comfort and user weight
- Removable and replaceable four-way stretch polurethane cover, featuring a low moisture transmission rate, keeping the patient cool, dry and comfortable
- No-gap seating surface
- Water resistant to prevent liquid ingress
- Fire retardant fabric and foam
- Built in pressure adjustment pumps to adjust pressure for comfort or to refill air cells
- Completely adjustable for patient positioning, with ergonomic controls for caregiver
- Reclining seat back adjustment
- Seat tilt adjustment – tilt in space
- Leg rest tilt adjustment
- Head rest is height adjustable and removeable
- Multi-position fold away arm rests, wings and footrest pads provide ease of access and transfer
- Individual locking, non-marking castors
- Single rear directional lock castor

Included accessories:
Headrest, Footrest Pads

Optional accessory:
Meal tray

Notice: all products, product specifications and data are subject to change without notice to improve reliability, function, design or otherwise

info@pegasushealthcare.com.au
03 9596 6377

437 Nepean Highway,
Brighton East, 3187 Australia

www.pegasushealthcare.com.au

