

Patient Safety and Comfort

- 3 cell cycle alternation with a 15 minute cycle time provides regular relief of interface pressure
- Wide diameter tubing enables rapid air-cell deflation 4 times per hour which stimulates reactive hyperaemia in the deep tissues surrounding the vulnerable bony prominences
- Cell-in-cell technology improves patient comfort and safety by providing a permanently inflated base portion beneath the active alternating component of the cell
- Fully automatic adjustment to patient weight and body position
- Dual comfort settings available within both Static and Dynamic modes
- Dual pressure settings improve comfort for frail and underweight patients without compromising pressure relief
- Static mode provides a firm and stable surface for essential nursing procedures such as physiotherapy, wound dressing and changing bed linen
- Transport mode ensures the mattress remains inflated for 8 hours in the event of a mains power failure
- V-Lock system prevents accidental disconnection of the power cord

Additional Nursing Support

- Audible and visual alarms provide peace of mind for nursing staff
- Rapid Cardiopulmonary Resuscitation (CPR) deflation to suit individual patient comfort or clinical requirements
- Static mode defaults to alternating after 30 minutes
- For extra long and adjustable length hospital beds.

Maximises Infection Control

- Waterproof, vapour permeable, Dartex® covers provide microclimate control at the skin/mattress interface
- Soft, easy to clean cover contains an antimicrobial agent that can withstand high temperature washing
- Welded seams prevent the ingress of fluids into the mattress
- 360° zip ensures cover can be easily removed for infection control and laundering
- Generous waterfall flap protects the zip

Technical Specifications (Trio Extra Long)

Number of cells	20
Mattress sizes	2100 x 880 x 200mm
Max. patient weight	220 kg
Mattress material	Cells: Polyurethane Cover: Dartex® two-way stretch Polyurethane coated Nylon
Pump unit dimensions	400 x 160 x 250 mm
Mattress weight	9.2 kg
Pump weight	2.9kg
Classification	Class B Medical Device / Class II
Mattress - flammability rating	Complies with BS7175 ignition source 0, 1 and 5
Power unit - flammability rating / ingress protection	UL 94 V-0 / IP 42
Power supply / power consumption:	230 V~ ±10%, 50 / 60 Hz 20 W max

Australian Distributor:



Pegasus Healthcare

437 Nepean Highway
Brighton East
Victoria, 3187

Tel: 03) 9596 6377

info@pegasushealthcare.com.au

www.pegasushealthcare.com.au

www.roberlimited.com

Tel. +44 (0) 1246 245404
@ info@roberlimited.com
@RoberLtd



© 2016 AirFlex® is a registered trademark of Rober Limited



Rober Limited, Wheatbridge Road, Chesterfield, Derbyshire, S40 2AB, UK.



Rober Ltd reserves the right to make technical changes.