VICAIR® CENTRE RELIEF 02

OFFLOADING PRESSURE IN THE HIGH RISK AREA



For completely offloading of the rectal area

Adjustable non-contact area Stable positioning

The Vicair Centre Relief O2 is the best choice if you need to completely offload the rectum perineum and coccyx area. By simply removing SmartCells from the centre circle compartment you create a non-contact area. This important feature makes a great difference for persons with (for instance) localized deep pressure ulcers, rectum or tailbone problems. Or after perineal surgery.

100% machine washable and breathable, for optimal hygiene and perfect temperature and moisture regulation.

The Vicair Centre Relief O2 helps you to sit comfortably and safe, even in the most painful situations.

AAAA

Maintenance & Safety	
Durability	***
Reliability	***
Functionality when misused	
Low maintenance level	AAAA
Level of hygiene	***

Functionality & Application	
Comfort	AAAA
Low cushion weight	AAAA
Available size range	
Maximum user weight	AAAA
Low seat to floor height	
Auto adjustments	

USED WITH

CUSHION SIZES

Rectum problems, tailbone problems, localized deep pressure ulcers, after perineal surgery or difficult childbirth

TECHNICAL SPECIFI	CATIONS
Height	10 cm
Average weight	750 grams
Max. user weight	250 kg
Machine washable	60°C max.

	Star	nda	rd									
	On ı	equ	ues	t								
												,
Inch	Inch		15	16	17	18	19	20	21	22	23	
	СМ	35	38	40	43	45	48	50	53	55	58	,
14	35											
15	38											
16	40											
17	43											
18	45											
19	48											
20	50											
21	53											
22	55											
23	58											
24	60											